

Instructions for Use**RealLine HIV quantitative - Uni-Format**

KIT FOR THE EXTRACTION AND QUANTITATIVE DETECTION OF RNA FOR HIV (HUMAN IMMUNO-DEFICIENCY VIRUS) BY REAL TIME PCR









Attention!
Please read the information
about quantification
process carefully!

Research Use Only (RUO)

RealLine HIV quantitative (Uni format) incl. Extraction	VBD0199	48 Tests
valid from	December 2019	

RealLine HIV quantitative Uni-Format

Explanation of symbols used in labeling

	For research use only
	Batch code
REF	Catalogue number
	Content of number of tests
	Expiry date
	Temperature limitation
	Consult instructions for use
	Manufacturer
	Keep out of sunlight



BIORON Diagnostics GmbH

In den Rauhweiden 20
67354 Römerberg (Germany)
Phone +49 6232 298 44 0
Fax +49 6232 298 44 29
info@bioron.de

Trademarks:

Rotor-Gene® is a registered trademark of Qiagen Group, Germany.

BIORON Diagnostics GmbH
info@bioron.de - www.bioron.de

RealLine HIV quantitative Uni-Format

Table of Content:

1. INTENDED USE	4
2. KIT CONTENTS	5
3. PRINCIPLE OF THE METHOD	6
4. SPECIFICATIONS	7
5. PRODUCT USE LIMITATIONS	7
6. WARNINGS AND PRECAUTIONS	8
7. ADDITIONAL MATERIALS AND DEVICES REQUIRED BUT NOT SUPPLIED	8
8. SPECIMEN TRANSPORT AND STORAGE	9
9. REAGENT PREPARATION	10
10. PROCEDURE PROTOCOL	11
10.1. UltraSensitive procedure for RNA isolation: from 1 ml of serum (<i>plasma</i>)	11
10.2. Standard protocol for RNA isolation: from 100 µl of serum (<i>plasma</i>):	12
10.3. PCR	14
11. DATA ANALYSIS AND INTERPRETATION	15
12. STORAGE AND TRANSPORTATION	16
ANNEX I: CALCULATING THE CONCENTRATION OF SPECIFIC VIRAL NA	17
ANNEX II: Programming the device and analysis of results using Rotor-Gene cyclers:	18
ANNEX III PROGRAMMING OF DT-96 CYCLER AND REALLINE CYCLER :	19

RealLine HIV quantitative Uni-Format

EXTRACTION AND QUANTITATIVE DETECTION OF HIV RNA BY REAL TIME PCR

Research Use Only

1. INTENDED USE

The assay kit **RealLine HIV quantitative** is intended for detection of detection and quantitative determination of human immunodeficiency virus (HIV) RNA in blood serum (plasma) using a method based on a reverse transcription of viral RNA followed by an amplification of cDNA in polymerase chain reaction (RT-PCR) with real-time hybridisation fluorescence detection.

The kit is intended for use with the recording thermocyclers of both - rotor type: Rotor-Gene 3000, Rotor-Gene 6000 (Corbett Research, Australia), Rotor-Gene Q (Qiagen, Germany), and microplate type: iQ iCycler, iQ5 iCycler, CFX96 (Bio-Rad, USA), DT-96 (DNA-Technology, Russia), RealLine Cyclers (BIORON Diagnostics GmbH) or analogues thereof.

The assay kit contains reagents sufficient for 4 x 12-test runs, which may be performed separately or simultaneously.

The assay kit can be used with either of two specimen preparation procedures, the Standard Procedure or the UltraSensitive procedure.

- In the Standard specimen preparation procedure, HIV RNA is isolated from **100 µl** of serum (plasma).
- In the UltraSensitive specimen preparation procedure, HIV viral particles in serum (plasma) are concentrated by Concentrating Solution of **1 ml** of serum (plasma).

RealLine HIV quantitative Uni-Format

2. KIT CONTENTS

Specimen Preparation Reagents:	
Concentrating Solution	4 vials, 14 ml each
Lysis Reagent № 1	4 vials, 4 ml each
Lysis Reagent № 2	4 vials, 7 ml each
Sorbent (suspension of magnetic particles)	1 vial, 1 ml
Solution for NA Precipitation	4 vials, 12 ml each
Wash Solution № 1	4 vials, 8 ml each
Wash Solution № 2	4 vials, 5 ml each
Specimen Diluent	4 vials, 3 ml each
Control samples:	
Recovery Solution for Control samples (RSC)	2 vials, 4 ml each
Positive Control (PC) , lyophilized	2 vials
Negative Control (NC) (Negative serum from human)	2 vials, 12 ml each
Internal Control (IC) , lyophilized	2 vials
CS1 and CS2: Samples for calibration, lyophilized - are used when the adequacy of analytical system has to be checked, see the Attachment 2	1 vial of each
Amplification reagents	
Ready Master Mix for reverse transcription and PCR, lyophilized (RMM)	48 test tubes
Additional	
The kit also contains plastic caps for vials with control and calibration samples	6
Passport with the specific amounts for the controls in the upper lid of the box	

RealLine HIV quantitative Uni-Format

3. PRINCIPLE OF THE METHOD

The kit contains all the required components to perform the whole range of procedures to concentrate and isolate human immunodeficiency virus RNA from blood serum (*plasma*), perform PCR with hybridisation fluorescent detection in real time, and calculate the amount of human immunodeficiency virus RNA taking into account the calibration and control samples.

The testing methodology is based on recording the amplification of a selected specific fragment of cDNA (*obtained by reverse transcription from HIV RNA*) that involves repeating cycles of: temperature denaturation, annealing of primers with complementary sequences, extension of polynucleotide sequences from these primers with Taq polymerase.

The employed detection method is based upon measuring fluorescence signals in each PCR cycle. The enhancement of fluorescence signals occurs due to the use of a hybridisation DNA probe that is specific for this cDNA, which binds to one of the DNA strands during the reaction, thus ensuring additional specificity of the method. The DNA-probe contains fluorescent dye at the 5' end and fluorescent quencher, which considerably decreases the fluorescence intensity, at the 3' end. In the course of complementary chain polymerase synthesis, due to 5'-3' nuclease activity of Taq DNA polymerase, the probe splits at the 5' end, and disintegration of dye and quencher takes place, which leads to increasing fluorescence signal as the reaction product is accumulated. As this takes place, the measured fluorescence intensity depends on the amount of generated specific amplicons, and the dynamics of fluorescence increase is determined by the initial amount of cDNA in the sample.

Quantitative evaluation of the efficiency of RNA isolation from the samples is ensured by the isolation of human immunodeficiency virus RNA from the blood serum (plasma) samples along with the preliminarily added internal control sample (IC), which allows controlling the efficiency of RNA isolation from the test samples and taking the same into account in calculation of the HIV RNA concentration.

RealLine HIV quantitative Uni-Format

4. SPECIFICATIONS

- 4.1. **Specificity** of human immunodeficiency virus RNA (*using negative sera of the manufacturer's Standard Reference Panel*) is 100%.
- 4.2. The minimum reliably detected HIV RNA concentration is 20 IU/ml.
- 4.3. The linear range is from 100 to 1×10^8 IU/ml.
- 4.4. The amplification efficiency is (100 ± 10) %.
- 4.5. The variation coefficient calculated for logarithmic values of human immunodeficiency virus RNA concentration is no more than 10%.
- 4.6. The linearity test, which characterises coincidence between the calculated value and the prescribed value (calculated taking into account the dilution factor) of HIV RNA concentration logarithm in the sample prepared from a sample with the known HIV RNA concentration, is (100 ± 10) %.













5. PRODUCT USE LIMITATIONS

For Research Use Only.

Not for use in diagnostic procedures. No claim or representation is intended to provide information for the diagnosis, prevention, or treatment of a disease.

RealLine HIV quantitative Uni-Format

6. WARNINGS AND PRECAUTIONS

-  For in vitro use only.
-  The kits must be used by skilled personnel only.
-  When handling the kit, follow the national safety requirements for working with pathogens.
-  To prevent contamination, the stages of DNA isolation and PCR test run must be spatially separated.
-  Avoid microbial and nuclease contamination of reagents when removing aliquots from reagent vials.
-  Wear protective disposable gloves, laboratory coats and eye protection when handling specimens and kit reagents.
-  Every workplace must be provided with its own set of variable-volume pipettes, necessary auxiliary materials and equipment. It is prohibited to relocate them to other workplaces.
-  The use of sterile disposable pipette tips is recommended.
-  Never use the same tips for different samples.
-  Do not pool reagents from different lots or from different vials of the same lot.
-  Dispose unused reagents and waste in accordance with country, federal, state and local regulations.
-  Do not use the kit after the expiration date at the side label of the kit.

7. ADDITIONAL MATERIALS AND DEVICES REQUIRED BUT NOT SUPPLIED

- Real time PCR device RNA Extraction Kit: RealLine Extraction 100 or RealLine Extraction 1000
- RealLine Internal Control (VBC8881) and Negative Control, if Extraction Kit from other supplier is used
- Tubes or plates and PCR foils, if the delivered format is not matching with the used system.
- Microcentrifuge (min RCF 13000 rpm);
- Vortex mixer;
- Thermo Shaker
- Disposable gloves, powder-free;
- Pipettes (capacity 10-100 µl, 100-1000 µl) with filters (aerosol barriers);
- Disposable DNase/RNase-free tips with filters
- Displacement tips
- 2.0 ml polypropylene tubes, sterile, non-siliconised
- Vacuum aspirator;
- 2.0 ml and 0.2 ml microtube racks;
- Magnetic Rack for nucleic acids isolation (e.g. BIORON Diagnostics GmbH).

RealLine HIV quantitative Uni-Format

8. SPECIMEN TRANSPORT AND STORAGE

Preparation of blood serum (plasma)

Attention! *Samples anticoagulated with heparin are unsuitable for the test.*

Transfer serum (plasma) to a sterile polypropylene tube within 6 hours after blood sampling.

Transportation and storage of samples:

- At (18 – 25) °C – for no more than 2 hours;
- At (2 - 8) °C – for no more than 24 hours
- At minus 18 °C and below – up to 2 weeks

Only one freeze-thaw cycle is allowed!

Prior to use, centrifuge samples of serum (plasma) at (18 – 26) °C, 13,000 rpm for 5 min.

RealLine HIV quantitative Uni-Format

9. REAGENT PREPARATION

9.1. Prior to work, take the assay kit from the refrigerator, open the package, and keep the kit components at room temperature (18 - 25) °C for 30 minutes.

9.2. Open the vial with **PC** by removing plastic cap and rubber stopper. Put the removed caps and stoppers into the container with disinfecting solution.

Add 1 ml of RSC, tightly close the vial with new plastic cap provided with the kit. Carefully mix, keep at room temperature for 15 minutes, and then thoroughly mix again.

After diluting store at (2–8)°C for no longer than 1 month.

9.3. Open the vials with **CS1** and **CS2** by removing plastic caps and rubber stoppers. Put the removed caps and stoppers into the container with disinfecting solution.

Add **1 ml of RSC**, tightly close the vials with new plastic caps provided with the kit. Carefully mix, keep at room temperature for 15 minutes, and then thoroughly mix again

Store diluted CS1 and CS2 store at (2 – 8) °C for no more than 1 month.

9.4. Open the vial with **IC** by removing plastic cap and rubber stopper. Put the removed caps and stoppers into the container with disinfecting solution.

Add **1 ml of RSC**, tightly close the vial with new plastic cap provided with the kit. Carefully mix, keep at room temperature for 15 minutes, and then thoroughly mix again.

Store diluted IC at (2 – 8) °C for no more than 1 month.

Negative Control (NC) sample is ready to use.

Attention! *Once opened, NC should be stored at (2 - 8) °C and used within 1 month.*

9.5. Prior to use, heat **Lysis reagent № 1 and № 2** at (50 - 60) °C and mix thoroughly to dissolve the precipitated material. Vortex Sorbent to a condition of homogeneous suspension. Add **140 µl** of Sorbent suspension into a vial with Lysis Reagent № 2. Mix carefully.

Attention! *Lysis Reagent No. 2 is not to be stored after opening and adding Sorbent.*

For the quantitative determination, the RNA extraction is conducted from 100 µl of serum (plasma) using this kit. The Internal Control must be used with the Extraction, if another RNA Extraction kit is used, please add **RealLine Internal Control (REF VBC8881)** to the extraction procedure

Each set of samples must contain 3 PC and 1 NC sample.

The elution volume have to be 200 µl.

Please use for the NC samples the Negative Control Sample provided with the kit and do not use water or buffer.

RealLine HIV quantitative Uni-Format

10. PROCEDURE PROTOCOL

It is recommended to use three replicas of Positive Control sample and one Negative Control sample in each test run.

10.1. UltraSensitive procedure for RNA isolation: from 1 ml of serum (*plasma*)

10.1.1 When choosing the number of test samples for RNA isolation ensure that each group of isolated samples and each individual PCR run contain 3 CS1, 3 CS2, 1 NC, and 1 PC.

Attention! CS1 and CS2 undergo the RNA extraction procedure and RT-PCR during the first reaction with the kits of this lot and cycler in use. Further RT-PCR test runs can be done using only two control points: NC and PC, the calibration plot can be imported from the first analysis.

10.1.2 Prepare and label an appropriate number of 2 ml tubes for the test samples, PC, and NC.

10.1.3 Pipette **30 µl of IC** to each tube.

10.1.4 For NC, add to the tube, marked NC, **1 ml of Negative Control NC** sample.

10.1.5 For PC, add to the tube, marked PC, **970 µl of NC**, and **30 µl of PC sample**

10.1.6 For CS1 and CS2, pipette **970 µl of NC** and **30 µl** of the corresponding **calibration sample**

10.1.7 Add **1 ml** of each specimen to the appropriately labeled tube.

10.1.8 For each specimen or control tube, add **1 ml of Concentrating Solution**.

10.1.9 Close the tubes with caps, thoroughly mix the contents by rolling over the tubes 5 times, and keep for 10–15 minutes at room temperature. Then spin at 3,000 rpm at room temperature for 5 min*

* The spinning parameters depend on the model of the centrifuge in use:

Eppendorf MiniSpin — 3,000 rpm;

Eppendorf 5415D — 2,000 rpm;

Eppendorf 5417C — 2,000 rpm; and

Heraeus BiofugePico — 2,000 rpm.

10.1.10 Using a new tip for each sample, carefully remove the supernatant without disturbing the pellet. The pellet should be clearly visible at this step.

10.1.11 Add **200 µl of Lysis reagent № 1** to each tube. Vortex vigorously for 10-15 seconds. Some insoluble material may remain. Leave tubes for 5 minutes at room temperature. Briefly spin to collect drops from the tube walls..

RealLine HIV quantitative Uni-Format

- 10.1.12** Add **500 µl of Lysis reagent No. 2** with sorbent to each tube (with Lysis reagent No. 1). Vortex the contents of the tubes for 10–15 seconds. Keep in thermal shaker at 56°C for 10 min at 1,300 rpm. Briefly spin to dump collect from the tube walls.
- 10.1.13** Add **750 µl of NA precipitator** into each tube with the test samples.
- 10.1.14** Vortex the contents of the tubes for 10–15 seconds. Keep for 3–5 minutes at room temperature (18–25) °C. Spin for 5 minutes at 13,000 rpm.
- 10.1.15** Without shaking up the pellet place the tubes in the magnetic stand. Using pipette (*or aspirator*) with individual tip, remove supernatant from each tube without touching the pellet.
- 10.1.16** Add **500 µl of Wash Solution No. 1** to the pellet in each tube. Vortex the contents of the tubes for 10–15 seconds. Spin at 13,000 rpm for 5 min.
- 10.1.17** Without shaking up the pellet place the tubes in the magnetic stand. Using pipette (*or aspirator*) with individual tip, remove supernatant from each tube without touching the pellet.
- 10.1.18** Add **300 µl of Wash Solution No. 2** to the pellet in each tube. Vortex the contents of the tubes for 10–15 seconds. Spin at 13,000 rpm for 5 min.
- 10.1.19** Without shaking up the pellet place the tubes in the magnetic stand. Using pipette (*or aspirator*) with individual tip, remove supernatant from each tube without touching the pellet.
- 10.1.20** Dry the pellets in open tubes at room temperature (18 – 25) °C for 2–3 min.
- 10.1.21** Add **200 µl of Specimen Diluent** to the pellet in each tube. Thoroughly resuspend the pellet by vortexing. Incubate in thermal shaker at 56°C for 10 min at 1,300 rpm. Centrifuge at 13,000 rpm for 1 min. Samples are now prepared for a PCR test run.

Attention! *Extracted RNA cannot be stored! Perform RNA extraction immediately before running the RT-PCR.*

10.2. Standard protocol for RNA isolation: from 100 µl of serum (plasma):

- 10.2.1** With small volume of test serum (*plasma*), the kit can be used to detect HIV RNA in **100 µl sample**. *Be aware that, in this case, detection of HIV RNA at concentrations lower than 200 IU/ml is not guaranteed.*
- 10.2.2** Pick out and sign the required number of **2 ml** tubes (*according to the number of test samples, including the required controls* — 3 CS1, 3 CS2, 1 NC, and 1 PC).
- 10.2.3** Pipette **30 µl of IC solution** into each tube.

RealLine HIV quantitative Uni-Format

- 10.2.4** For the negative control, pipette **100 µl of NC** into the tube (*labelled as NC*).
- 10.2.5** For the positive control, pipette **70 µl of NC** and **30 µl of PC** into the tube (*labelled as PC*).
- 10.2.6** Pipette **70 µl of NC** and **30 µl** of the corresponding **calibration sample** into the tubes labelled as **CS1** and **CS2**.
- 10.2.7** In other tubes, pipette **100 µl of the test blood serum or plasma samples** according to the tube labels. Use individual pipette tip for each sample.
- 10.2.8** Add **500 µl of Lysis reagent No. 2 with sorbent** to each tube. Vortex the contents of the tubes for 10–15 seconds. Keep in thermal shaker at 56°C for 10 min at 1,300 rpm. Briefly spin to dump drops from the tube walls.
- 10.2.9** Add **600 µl of NA precipitator** into each tube with the test samples.
- 10.2.10** Vortex the contents of the tubes for 10–15 seconds. Keep for 3–5 minutes at room temperature (18 – 25) °C. Spin for 5 minutes at room temperature at 13,000 rpm.
- 10.2.11** Without shaking up the pellet place the tubes in the magnetic stand. Using pipette (*or aspirator*) with individual tip, remove supernatant from each tube without touching the pellet.
- 10.2.12** Add **500 µl of Wash Solution No. 1** to the pellet in each tube. Vortex the contents of the tubes for 10–15 seconds. Spin at 13,000 rpm for 5 min.
- 10.2.13** Without shaking up the pellet place the tubes in the magnetic stand. Using pipette (*or aspirator*) with individual tip, remove supernatant from each tube without touching the pellet.
- 10.2.14** Add **300 µl of Wash Solution No. 2** to the pellet in each tube. Vortex the contents of the tubes for 10–15 seconds. Spin at 13,000 rpm for 5 min.
- 10.2.15** Without shaking up the pellet place the tubes in the magnetic stand. Using pipette (*or aspirator*) with individual tips, remove supernatant from each tube without touching the pellet.
- 10.2.16** Dry the pellets in open tubes at room temperature (18 – 25) °C for 2–3 minutes.
- 10.2.17** Add **200 µl of Specimen Diluent** to the pellet in each tube. Thoroughly resuspend the pellet by vortexing. Incubate in thermal shaker at 56°C for 10 min at 1,300 rpm. Spin at 13,000 rpm for 1 min. Samples are now prepared for a PCR test run.

Attention! *Extracted RNA cannot be stored! Perform RNA extraction immediately before running the RT-PCR.*

RealLine HIV quantitative Uni-Format

10.3. PCR

10.3.1 Place the tubes with prepared samples and controls to the magnetic Rack.

10.3.2 Prepare an appropriate number of reaction tubes with Ready Master Mix (RMM). Label each reaction tube for each specimen and control.

Attention! For blockcycler mark tubes on the lateral part, for Rotor-Gene mark on the cap.

10.3.3 Add **50 µl** of corresponding extracted RNA solution and control to the appropriately labeled reaction tube using a new RNase-free tip with aerosol barrier for each sample. **Do not grasp sorbent particles!**

10.3.4 Place reaction tubes into the thermal block of real time PCR device and program real time PCR device as follows:

For iQ/iQ5 iCycler, CFX96, DT-96, RealLine Cycler:

Step 1:	45°C	30 min	
Step 2:	94°C	1 min	
Step 3:	94°C	10 sec	50 Cycles
	60°C*	20 sec	
* measurement of fluorescent at 60 °C in FAM and ROX			

For Rotor-Gene 3000/6000/Q:

Step 1:	45°C	30 min	
Step 2:	94°C	1 min	
Step 3:	94°C	10 sec	50 Cycles
	60°C*	40 sec	
* measurement of fluorescent at 60 °C in FAM and ROX			

10.3.5 Select the amplification detection channels:

- collect data through the **FAM** channel (iQ5 iCycler, CFX96, DT-96, Rotor-Gene 3000) or the **Green** channel (Rotor-Gene 6000, Rotor-Gene Q) for the detection of amplification signal of **IC cDNA**;
-
- collect data through the **ROX** channel (iQ5 iCycler, CFX96, DT-96, Rotor-Gene 3000) or the **Orange** channel (Rotor-Gene 6000, Rotor-Gene Q) for the detection of amplification signal of **HIV cDNA**

10.3.6 Program the position of tubes with samples, calibration samples, positive and negative control samples according to the instruction manual to the used device.

10.3.7 Run the program and perform PCR with fluorescent detection in real time.

RealLine HIV quantitative Uni-Format

11. DATA ANALYSIS AND INTERPRETATION

Perform analysis, evaluation of results, and calculation of **HIV RNA** concentrations in the samples according to the instruction manual to the used device.

11.1. General conditions of analysis and evaluation of results

11.1.1. In **PC**, the following should be registered:

- an increasing signal of the specific amplification product for **HIV cDNA** (the **ROX (Orange)** channel) and determined the threshold cycle value, **Ct**, for **HIV**;
- an increasing amplification signal for **IC cDNA** in **FAM (Green)** channel and determined the threshold cycle value, **Ct**, for **IC**.

11.1.2. The results of an individual RT-PCR run are to be analysed and evaluated if, for **PC**, the **HIV Ct** value is within the range indicated in the insert label for the kits of this lot.

11.1.3. For **NC**, an increasing amplification signal for **IC DNA** should be registered and **IC Ct** should be determined with no significant increase in the signal for the specific amplification product of **HIV RNA**. No **ROX (Orange)** fluorescence increase should appear.

11.1.4. For each **sample**, the program should detect the increase in the amplification signal of **IC cDNA** in **FAM (Green)** channel and determine **IC Ct**. The result of the test sample is considered valid if **IC Ct in FAM (Green)** is **less than or equal to 40**.

11.2. Evaluation of results

11.2.1. Calculate **(IC Ct)_{av}** as an average value for **IC Ct** for all test samples (*including controls*). The **IC Ct** values differing by more than 2 from the **(IC Ct)_{av}** value should be discarded. After discarding, recalculate **(Ct IC)_{av}** for the remaining values.

11.2.2. The test sample is considered as **positive**, i.e. containing **HIV RNA**, if, for this sample, the **Ct** value in the **ROX (Orange)** channel is **less than or equal to 40**.

11.2.3. The test sample is considered as **negative** (*not containing HIV RNA*) if, for this sample, the **Ct** value in the **ROX (Orange)** channel is **above 40 or not determined**.

If, for such sample, the **IC Ct** value exceeds the **(IC Ct)_{av}** value by more than 2, the result for this sample is not to be analysed and is evaluated as **negative**. A repeated test for this sample starting from the isolation step is required.

In case of repetitive result take another blood sample and repeat the testing once again.

11.2.4. In case of contamination **all positive** results for this individual PCR run are **considered invalid**. Take measures to detect and eliminate the contamination source and repeat the testing

RealLine HIV quantitative Uni-Format

for all samples of this run, for which the positive result has been obtained. The samples of this run, for which the test has yielded negative results, should be considered as negative.

11.3. Quantitative analysis of results

For quantification, calculate HIV RNA concentration in the test samples according to **Annex I**.

Attention! If DNA was isolated from 100 µl of sample, multiply the result by 10.

11.3.1. If the calculated value of HIV RNA concentration in the sample is within the range of **100 IU/ml to 10⁸ IU/ml**, the result is considered as **positive**, with indication of the obtained HIV RNA concentration in the sample (*in IU/ml*).

11.3.2. If the calculated value of HIV RNA concentration in the sample is higher than 10⁸ IU/ml, the result is recommended to be interpreted as **positive**, with HIV RNA concentration “**more than 10⁸ IU/ml**”.

11.3.3. If the calculated value of HIV RNA concentration is less than 100 IU/ml, the result is recommended to be interpreted as **positive**, with HIV RNA concentration “**less than 100 IU/ml**”.

11.3.4. The test sample is considered **negative** (*not containing HIV RNA*) if, for this sample, if the **Ct** value in the **ROX (Orange)** channel is **above 40 or not determined**.

11.3.5. If the obtained concentration of PC is not within the range indicated in the insert label for the kit, repeat the tests for all samples from this run.

12. STORAGE AND TRANSPORTATION

- Store the assay kit at (2 - 8) °C in the manufacturer's packing.
- Transport at (2 – 8) °C. Transportation for up to 25 °C for up to 10 days is allowed.
- Do not freeze the kit!
- Do not pool reagents from different lots or from different vials of the same lot.
- Strictly follow the Instruction manual for reliable results.
- Do not use kits with damaged inner packages and get in contact with BIORON Diagnostics GmbH.

RealLine HIV quantitative Uni-Format

ANNEX I: CALCULATING THE CONCENTRATION OF SPECIFIC VIRAL NA

Calculation of HIV RNA concentration

Quantitative analysis of HIV RNA is performed using the formulas below.

- 1) Calculate coefficient β :

$$\beta = [\lg(C_{CS1}) - \lg(C_{CS2})] / [\Delta Ct_{CS2} - \Delta Ct_{CS1}],$$

where C_{CS1} and C_{CS2} are the concentrations of specific viral NA in the calibration samples (*indicated in the insert label for the kit*);

$$\Delta Ct_{CS1} = Ct_{CS1av} - IC\ Ct_{CS1av};$$

$$\Delta Ct_{CS2} = Ct_{CS2av} - IC\ Ct_{CS2av};$$

Ct_{CS1av} , Ct_{CS2av} are the average Ct values for **CS1** and **CS2** samples through the “**ROX**” channel;

$IC\ Ct_{CS1av}$, $IC\ Ct_{CS2av}$ are the average Ct values for **CS1** and **CS2** samples through the “**FAM**” channel.

- 2) Calculate the concentration of specific viral NA in the test sample (C_k):

$$C_k = C_{CS2} \times 10^{\beta(\Delta Ct_{CS2} - \Delta Ct_k)},$$

where k is the sample number;

$$\Delta Ct_k = Ct_k - IC\ Ct_k$$

where Ct_k and $IC\ Ct_k$ are the Ct values for the sample number k in the “**ROX**” and “**FAM**” channels, respectively.

If extraction and RT-PCR test run for **CS1** and **CS2** are not performed, use the β and ΔCt_{CS2} values obtained in the first test run with the kit of this lot for calculation.

For an easier calculation we can provide you with an Excel sheet. Ask us:

techsupport@bioron.de

Note: For precise calculation of viral NA concentration in analyzed samples, the validation of the analytical system should be done by comparison of **Positive Control PC** concentration considering Internal control concentrations. The PC concentration is specified in the passport of the assay kit. For this purpose three independent Positive Control samples are needed.

RealLine HIV quantitative Uni-Format

ANNEX II: Programming the device and analysis of results using Rotor-Gene cyclers:

Rotor-Gene 3000, Rotor-Gene 6000 (Corbett research, Australia), Rotor-Gene Q (Qiagen, Germany)

Hereinafter, detection channels and terms corresponding to different versions of devices and software are listed in the following order: Rotor-Gene 3000 (Rotor-Gene 6000, Rotor-Gene Q).

- Click **New** button.
- Select an **Advanced** template from the tab of the New Run wizard. Click **New** button.
- Select **36-Well Rotor** type, check that No Domed 0.2 ml Tubes are used. Click **Next** button.
- In the new window, determine Reaction volume as **50 µl**. Click **Next** button.
- The temperature profile of real time PCR should be set. Click **Edit Profile** button.

• Step 1:	• 50°C	• 2min	•
• Step 2:	• 95°C	• 2min	•
• Step 3:	• 94°C	• 10 sec	• 50 cycles
•	• 60°C*	• 40 sec	

 - * Measure the fluorescence at 60°C
- Then temperature profile is set, click **OK** button.
- In the **New Run Wizard** window click **Calibrate (Gain optimization)** button. The window **Auto Gain Calibration Setup** opens. In the line **Channel Settings** choose **ROX (Orange)**, click **Add**. Set **Tube Position 1**, **Min Reading 5**, **Max Reading 10**, click **OK**. In the line **Channel Settings** choose **FAM (Green)**, click **Add**. Set **Tube Position 1**, **Min Reading 5**, **Max Reading 10**, click **OK**.
- Tick off **Perform Calibration Before 1st Acquisition**. Click **Close** button.
- Click **Next** button, start the amplification process by clicking **Start Run** button.
- Save a file in the Rotor-Gene/templates folder, named RealLine with *.ret extension. In subsequent work RealLine template would be presented in New run wizard.
- Save reaction result file with Rotor-Gene Run File *.rex extension.
- Record the positions of the controls and specimens according to the instruction manual of the operating device. Click **Start run** button.

Results for Internal Control cDNA amplification

- Click **Analysis** button, choose **Quantitation** from the list, choose **Cycling A. FAM**, click **Show** button.
- Click **OK** button, and cancel automatic **Threshold** determination.
- Click **Linear scale** button.
- In the **Quantitation analysis** menu buttons **Dynamic tube** and **Slope Correct** should be pressed.
- Click **More Settings (Outlier removal)** button, determine **NTC threshold** value as **5%**.
- In the column **CT Calculation (right part of the window)** determine **Threshold** value as **0.04**.
- In the result table (**Quant. Results window**) **Ct** will be displayed.

RealLine HIV quantitative Uni-Format

Results for HIV cDNA amplification

- Click **Analysis** button, choose **Quantitation** from the list, choose **Cycling A. ROX** click, **Show** button.
- Click **OK** button, and cancel automatic **Threshold** determination.
- Click **Linear scale** button.
- In the **Quantitation analysis** menu buttons **Dynamic tube** and **Slope Correct** should be pressed.
- Click **More Settings (Outlier removal)** button, determine **NTC threshold** value as **5%**.
- In the column **CT Calculation** (*right part of the window*) determine **Threshold** value as **0.04**.
- In the result table (**Quant. Results window**) **Ct** will be displayed.

ANNEX III PROGRAMMING OF DT-96 CYCLER AND REALLINE CYCLER :

The measurement exposure must be adjusted. Choose the **Operation with the device** mode in the **Settings** menu, select the item **Measurement exposition:**

FAM to **500**;

HEX and **ROX** to **1000**;

Cy5 to **500**.

Confirm that the current exposure value is saved by pressing **YES**.

Attention! *The specified exposure values are applicable only for RealLine kits and, if necessary, should be changed for other purposes.*

RealLine HIV quantitative Uni-Format

